




Urologie
Urology
Urologie
Urología
Urologia

Prostata Haken, -Faßzange
Prostatic Retractors, Grasping Forceps
Ecarteurs à prostate, Pince à saisir la prostate
Separadores, Pinza para coger prostata
Divaricatori, Pinza per afferrare la prostata

ASH INSTRUMENTS

SURGICAL INSTRUMENTS



YOUNG
66101.21
21,5 cm




YOUNG
66102.22
21,5 cm



YOUNG
66103.22
20,5 cm



YOUNG
66105.21
22 cm



YOUNG
66106.20
20 cm

Kapsel­fass­zangen, Faden­föhrungs­zange
 Capsule Forceps, Ligature Carrying Forceps
 Pinc­es pour saisir les capsu­les. Pince porte-fil
 Pin­zas capsu­lares, Pin­ze porta ligaduras
 Pin­ze per afferrare capsu­le, Pin­za porta-filo



MILLIN
66111.23
 24,5 cm



MILLIN
66115.24
 24,5 cm




MILLIN
66117.23
 22 cm



MILLIN
66123.23
 22 cm

Fass­zange für vordere Lappen
 Lobe Holding Forceps
 Pince pour saisir les lom­beaux
 Pin­za para lóbulos anteriores
 Pin­za per lobuli anteriori




MILLIN
66125.24
 24,5 cm

Faden­föhrungs­zange
 Ligature Carrying Forceps
 Pince porte-fil
 Pin­ze porta ligaduras
 Pin­za porta-filo

Blasenwundsperrer, Blasenspatel
Bladder Retractors, Bladder Spatula
Ecarteurs vésicaux, Spatule vésicale
Separadores vesicales, Espátula vesical
Divaricatore vesicale, Spatole vescicali

ASH INSTRUMENTS

SURGICAL INSTRUMENTS



35 mm

LEGUEU
66147.27
27 cm



67 mm

KOCHER
66149.27
26 cm



JUDD-MASSON
66150.00



MILLIN

66151.50

MILLIN

66151.00

komplett mit je 2 Seitenvälven
57 x 25 mm und 80 x 25 mm,
sowie 1 Mittelvalve 120 x 45 mm
ohne Kaltlichtaufsatz

complete with 2 pairs of lateral
blades 57 x 25 mm and 80 x 25 mm
and 1 central blade 120 x 45 mm
without F.O. Carrier

complet, avec 2 paires de valves laterales,
57 x 25 mm et 80 x 25 mm,
ainsi qu'une valve centrale 120 x 45 mm
sans support lumiere froide

completo, con 2 pares de valvas laterales
57 x 25 mm y 80 x 25 mm
y una valva central 120 x 45 mm
sin portador de luz fria

completo, con 2 paia di valve laterali
57 x 25 mm e 80 x 25 mm
in piú 1 valva centrale 120 x 45 mm
senza supporto per luce fredda

MILLIN

66151.50

komplett mit je 2 Seitenvälven
57 x 25 mm und 80 x 25 mm,
sowie 1 Mittelvalve 120 x 45 mm
mit Kaltlichtaufsatz

complete with 2 pairs of lateral
blades 57 x 25 mm and 80 x 25 mm
and 1 central blade 120 x 45 mm
with F.O. Carrier

complet, avec 2 paires de valves laterales,
57 x 25 mm et 80 x 25 mm,
ainsi qu'une valve centrale 120 x 45 mm
avec support lumiere froide

completo, con 2 pares de valvas laterales
57 x 25 mm y 80 x 25 mm
y una valva central 120 x 45 mm
con portador de luz fria

completo, con 2 paia di valve laterali
57 x 25 mm e 80 x 25 mm
in piú 1 valva centrale 120 x 45 mm
con supporto per luce fredda

MILLIN

66151.51

Kaltlichtaufsatz
F.O. Light Carrier
Supporto lumiere froide
Portador de luz fria
Supporto per luce fredda



THOMSON-WALKER
66160.00

komplett mit
2 Paar Seitenvalven,
1 Mittelvalve und
1 Retraktor

complete with
2 pairs of lateral blades
1 central blade and
1 retractor

complet avec
2 paires de valves laterales
1 valve centrale
1 écarteur

completo con
2 pares de valvas laterales
1 valva central
1 separador

completo con
2 paia di valve laterali
1 valva centrale e
1 retrattore



LEGUEU
66163.00

Blasenhalsperrer
Bladder Neck Spreader
Ecarteur pour le col de la vessie
Dilatador para el cuello vesical
Divaricatore del collo vesciale

ASH INSTRUMENTS

SURGICAL INSTRUMENTS



MILLIN
66165.00
29 cm




MATHIEU
66168.20
Schaftlänge 20 cm
length of shaft 20 cm
longueur du tige 20 cm
longitud del vástago 20 cm
lunghezza del gambo 20 cm

Fremdkörperzangen
Foreign Body Forceps
Pinceps à corps étrangers
Pinzas para cuerpos extranjerios
Pinze per corpi estranei



MATHIEU
66168.28
Schaftlänge 28 cm
length of shaft 28 cm
longueur du tige 28 cm
longitud del vástago 28 cm
lunghezza del gambo 28 cm



OTIS
66169.00

komplett mit Holzetui
mit Dilator, 2 Messern und 3 Ansätzen

complete in wooden case
with dilator, 2 knives and 3 tips

complet, en boîte métal
avec dilateur, 2 lames et 3 embouts

completo an caja de madera
con dilatador, 2 cuchillos y 3 puntas

completo, in contenitore di legno
con dilatatore, 2 bisturi e 3 punte



OTIS
66169.01


Ersatmesser
Spare Blade
Couteau de rechange
Cuchillo de repuesto
Lama die ricambio

Circumcisionsgerät
 Circumcision Instrument
 Instrument de circoncision
 Instrumento para circuncisión
 Strumento per circuncisione



66171.08	Ø 8 mm
66171.11	Ø 11 mm
66171.13	Ø 13 mm
66171.16	Ø 16 mm
66171.21	Ø 21 mm
66171.26	Ø 26 mm
66171.32	Ø 32 mm

Urethrographieinstrument
 Urethrography Instrument
 Instrument pour urétrographie
 Instrumento para uretrografia
 Strumento per uretrografia



66172.01

KNUTSON-KOCH
66172.01

Konus
 Cone
 Cône
 Cono
 Cono

KNUTSON
66172.00

Beschneidungsinstrument
 Circumcision Instrument
 Instrument de circoncision
 Instrumento para circuncisión
 Strumento per circuncisione



MOGAN
66174.01
 für Erwachsene
 for adults
 pour adultes
 para adultos
 per adulti

MOGAN
66174.02
 für Kinder
 for children
 pour enfants
 para niños
 per bambini



WINKELMANN
66174.12 Ø 12 mm
66174.16 Ø 16 mm
66174.20 Ø 20 mm

Penisklemmen
Meatus Clamps
Compresseurs de l'urètre
Compressores uretrales
Serrapene



STOCKMANN
66175.08
8 cm



STRAUSS
66177.11
11 cm

Klemme zum Fixieren des Samenstranges
Fixation clamp for the spermatic cord
Pince pour fixer le cordon
Pinza para fijar el cordón espermático
Pinza per fissare il tronco spermatico

Zange zum Erfassen der Glans und der Vorhaut
Forceps for grasping the glans and the foreskin




66179.12
12,5 cm



66179.15
14,5 cm

Sekretlöffel, Katheterspanner
Secretions Scoops, Catheter Introducer
Cuiller à sécrétions, Introducteur de sondes
Cucharillas p. secreciones, Introducer de catéres
Cucchia per secrezioni, mandrini per sonde uretrali



GUTTMANN
66180.26
26,5 cm

DITTEL
66182.46
46 cm

GUYON
66183.46
46 cm

GUYON
66184.39
39 cm


GUYON
66185.35
35 cm

Gleitsonde, Harnröhrenbougies
Filiform Bougie, Urethral Sounds
Bougie conductrice, Bougies urétrales
Sonda filiforme, Bujias uretrales
Candeletta filiforme, Candelette dilatatrici

ASH INSTRUMENTS

SURGICAL INSTRUMENTS

mit Schraubgewinde-Ansatz
with screw thread
avec about fileté
con rosca de tornillo
con attacco filettato



FILIFORM

66199.03	3 Ch.
66199.04	4 Ch.
66199.05	5 Ch.
66199.06	6 Ch.


Satz mit 21 Bougies, Charr. 10 bis 30
mit 4 Gleitsonden
im Metalletui 42 x 17,5 x 4 cm

set with 21 bougies, Charr. 10 to 30
with 4 filiform bougies
in metal case 42 x 17,5 x 4 cm

jeu avec 21 bougies, Charr. 10 jusqu' à 30
avec 4 bougies conductrices
en boîte métal 42 x 17,5 x 4 cm

juego de 21 bujias, Charr. 10 a 30
con 4 bujias filiformes
en caja metálica 42 x 17,5 x 4 cm

set di 21 dilatatori, Charr. 10 fino a 30
con 4 candelette dilatatrici
in contenitore metallico 42 x 17,5 x 4 cm



GUYON
66200.00
compl.



VAN BUREN
66201.00
Set compl.
8-36 Charr.

van BUREN	mm Ø	Charrière (French)
66201.06	2	6
66201.08	2 ² / ₃	8
66201.10	3 ¹ / ₃	10
66201.12	4	12
66201.14	4 ² / ₃	14
66201.16	5 ¹ / ₃	16
66201.18	6	18
66201.20	6 ² / ₃	20
66201.22	7 ¹ / ₃	22
66201.24	8	24
66201.26	8 ² / ₃	26
66201.28	9 ¹ / ₃	28
66201.30	10	30
66201.32	10 ² / ₃	32
66201.34	11 ¹ / ₃	34
66201.36	12	36
66201.38	12 ² / ₃	38

mit Schraubgewinde-Ansatz
with screw thread
avec about fileté
con rosca de tornillo
con attacco filettato



GUYON
66203.00
Set compl.
10-30 Charr.


GUYON	mm Ø	Charrière (French)
66203.10	3 ¹ / ₃	10
66203.11	3 ² / ₃	11
66203.12	4	12
66203.13	4 ¹ / ₃	13
66203.14	4 ² / ₃	14
66203.15	5	15
66203.16	5 ¹ / ₃	16
66.203.17	5 ² / ₃	17
66203.18	6	18
66203.19	6 ¹ / ₃	19
66203.20	6 ² / ₃	20
66203.21	7	21
66203.22	7 ¹ / ₃	22
66203.23	7 ² / ₃	23
66203.24	8	24
66203.25	8 ¹ / ₃	25
66203.26	8 ² / ₃	26
66203.27	9	27
66203.28	9 ¹ / ₃	28
66203.29	9 ² / ₃	29
66203.30	10	30



GUYON

GUYON	mm Ø	Charrière (French)
66205.08	2 ² / ₃	8
66205.09	3	9
66205.10	3 ¹ / ₃	10
66205.11	3 ² / ₃	11
66205.12	4	12
66205.13	4 ¹ / ₃	13
66205.14	4 ² / ₃	14
66205.15	5	15
66205.16	5 ¹ / ₃	16
66205.17	5 ² / ₃	17
66205.18	6	18
66205.19	6 ¹ / ₃	19
66205.20	6 ² / ₃	20
62205.21	7	21
66205.22	7 ¹ / ₃	22
66205.23	7 ² / ₃	23
66205.24	8	24
66205.25	8 ¹ / ₃	25
66205.26	8 ² / ₃	26
66205.27	9	27
66205.28	9 ¹ / ₃	28
66205.29	9 ² / ₃	29
66205.30	10	30
66205.31	10 ¹ / ₃	31
66205.32	10 ² / ₃	32
66205.33	11	33
66205.34	11 ¹ / ₃	34
66205.35	11 ² / ₃	35
66205.36	12	36
66205.37	12 ¹ / ₃	37
66205.38	12 ² / ₃	38
66205.39	13	39
66205.40	13 ¹ / ₃	40
66205.50	16 ² / ₃	50

für Männer
 for men
 pour hommes
 para hombres
 per uomini



	mm Ø	Charrière (French)
66211.08	2 ² / ₃	8
66211.10	3 ¹ / ₃	10
66211.12	4	12
66211.14	4 ² / ₃	14
66211.16	5 ¹ / ₃	16
66211.18	6	18
66211.20	6 ² / ₃	20
66211.22	7 ¹ / ₃	22
66211.24	8	24

für Frauen
 for women
 pour femmes
 para mujeres
 per donne



	mm Ø	Charrière (French)
66213.08	2 ² / ₃	8
66213.10	3 ¹ / ₃	10
66213.12	4	12
66213.14	4 ² / ₃	14
66213.16	5 ¹ / ₃	16
66213.18	6	18
66213.20	6 ² / ₃	20
66213.22	7 ¹ / ₃	22
66213.24	8	24



COXETER	mm Ø	Charrière (French)
66219.08	2 ² / ₃	8
66219.10	3 ¹ / ₃	10
6219.12	4	12
66219.14	4 ² / ₃	14
66219.16	5 ¹ / ₃	16
66219.18	6	18
66219.20	6 ² / ₃	20
66219.22	7 ¹ / ₃	22
66219.24	8	24



WALTHER	mm Ø	Charrière (French)
66221.12	4	12
66221.14	4 ² / ₃	14
66221.16	5 ¹ / ₃	16
66221.18	6	18
66221.20	6 ² / ₃	20
66221.22	7 ¹ / ₃	22
66221.24	8	24
66221.26	8 ² / ₃	26
66221.28	9 ¹ / ₃	28
66221.30	10	30
66221.32	10 ² / ₃	32
66221.34	10 ¹ / ₃	34
66221.36	12	36
66221.38	12 ² / ₃	38

WALTHER
66221.00
 34,5 cm

Satz komplett
 complete set
 jeu complet
 juego completo
 set completo
66212.12 - 66212.38



DITTEL	mm Ø	Charrière (French)
66225.08	2 ² / ₃	8
66225.09	3	9
66225.10	3 ¹ / ₃	10
66225.11	3 ² / ₃	11
66225.12	4	12
66225.13	4 ¹ / ₃	13
66225.14	4 ² / ₃	14
66225.15	5	15
66225.16	5 ¹ / ₃	16
66225.17	5 ² / ₃	17
66225.18	6	18
66225.19	6 ¹ / ₃	19
66225.20	6 ² / ₃	20
66225.21	7	21
66225.22	7 ¹ / ₃	22
66225.23	7 ² / ₃	23
66225.24	8	24
66225.25	8 ¹ / ₃	25
66225.26	8 ² / ₃	26
66225.27	9	27
66225.28	9 ¹ / ₃	28
66225.29	9 ² / ₃	29
66225.30	10	30
66225.32	10 ² / ₃	32

DITTEL
66225.00

Satz komplett (23)
 complet set (23)
 jeu complet (23)
 juego completo (23)
 set completo (23)
66225.08 - 66225.30

DITTEL
66225.98

Satz komplett (6)
 complet set (6)
 jeu complet (6)
 juego completo (6)
 set completo (6)
66225.08/12/16/20/24/30



DITTEL
34,5 cm

DITTEL	mm Ø	Charrière (French)
66227.08	2 ² / ₃	8
66227.09	3	9
66227.10	3 ¹ / ₃	10
66227.11	3 ² / ₃	11
66227.12	4	12
66227.13	4 ¹ / ₃	13
66227.14	4 ² / ₃	14
66227.15	5	15
66227.16	5 ¹ / ₃	16
66227.17	5 ² / ₃	17
66227.18	6	18
66227.19	6 ¹ / ₃	19
66227.20	6 ² / ₃	20
66227.21	7	21
66227.22	7 ¹ / ₃	22
66227.23	7 ² / ₃	23
66227.24	8	24
66227.25	8 ¹ / ₃	25
66227.26	8 ² / ₃	26
66227.27	9	27
66227.28	9 ¹ / ₃	28
66227.29	9 ² / ₃	29
66227.30	10	30

